

Interleukin-1 receptor antagonist protein

Catalog No: #AP84787



Package Size: #AP84787-1 50ug #AP84787-2 100ug #AP84787-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor antagonist protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	IL-1RN,IL-1ra,IRAP,IL1 inhibitor
Accession No.	Uniprot ID:Q9BEH0
Uniprot	Q9BEH0
GeneID	403660;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL1RN

Note: This product is for in vitro research use only