

Keratin, type II cytoskeletal 1

Catalog No: #AP84788



Package Size: #AP84788-1 50ug #AP84788-2 100ug #AP84788-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type II cytoskeletal 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cytokeratin-1,CK-1,Epithelial keratin-1,Keratin-1,K1,Type-II keratin Kb1,KER1
Accession No.	Uniprot ID:Q6EIY9
Uniprot	Q6EIY9
GeneID	444857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT1

Note: This product is for in vitro research use only