

Keratin, type II cytoskeletal 6A

Catalog No: #AP84789



Package Size: #AP84789-1 50ug #AP84789-2 100ug #AP84789-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type II cytoskeletal 6A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin-6A,CK-6A,Keratin-6A,K6A,Type-II keratin Kb6
Accession No.	Uniprot ID:Q4FZU2
Uniprot	Q4FZU2
GeneID	683313;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt6a

Note: This product is for in vitro research use only