

cAMP-dependent protein kinase catalytic subunit alpha

Catalog No: #AP84792

Package Size: #AP84792-1 50ug #AP84792-2 100ug #AP84792-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	cAMP-dependent protein kinase catalytic subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PKA C-alpha,PKACA
Accession No.	Uniprot ID:P17612
Uniprot	P17612
GeneID	5566;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKACA

Note: This product is for in vitro research use only