

Arachidonate 5-lipoxygenase

Catalog No: #AP84794



Package Size: #AP84794-1 50ug #AP84794-2 100ug #AP84794-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Arachidonate 5-lipoxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	5-LO,LOG5
Accession No.	Uniprot ID:P09917
Uniprot	P09917
GeneID	240;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALOX5

Note: This product is for in vitro research use only