

## Integrin alpha-4

Catalog No: #AP84796



Package Size: #AP84796-1 50ug #AP84796-2 100ug #AP84796-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin alpha-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD49 antigen-like family member D, Integrin alpha-IV, Lymphocyte Peyer patch adhesion molecules subunit alpha, LPAM subunit alpha, VLA-4 subunit alpha, CD49d
Accession No.	Uniprot ID:Q00651
Uniprot	Q00651
GeneID	16401;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Itga4

Note: This product is for in vitro research use only