

Integrin beta-2

Catalog No: #AP84798



Package Size: #AP84798-1 50ug #AP84798-2 100ug #AP84798-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Complement receptor C3 subunit beta, CD18, CD18, MFI7
Accession No.	Uniprot ID:P05107
Uniprot	P05107
GeneID	3689;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ITGB2

Note: This product is for in vitro research use only