

## Inter-alpha-trypsin inhibitor heavy chain H1

Catalog No: #AP84799



Package Size: #AP84799-1 50ug #AP84799-2 100ug #AP84799-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inter-alpha-trypsin inhibitor heavy chain H1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ITI heavy chain H1
Accession No.	Uniprot ID:Q61702
Uniprot	Q61702
GeneID	16424;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Itih1

Note: This product is for in vitro research use only