

NF-kappa-B inhibitor beta

Catalog No: #AP84800



Package Size: #AP84800-1 50ug #AP84800-2 100ug #AP84800-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NF-kappa-B inhibitor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NF-kappa-BIB,I-kappa-B-beta,IkB-B,Thyroid receptor-interacting protein 9,TR-interacting protein 9,IKBB,TRIP9
Accession No.	Uniprot ID:Q15653
Uniprot	Q15653
GeneID	4793;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFKBIB

Note: This product is for in vitro research use only