

Submandibular glandular kallikrein-9

Catalog No: #AP84804



Package Size: #AP84804-1 50ug #AP84804-2 100ug #AP84804-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Submandibular glandular kallikrein-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rGK-9, KLK-S3, S3 kallikrein, Submandibular enzymatic vasoconstrictor, SEV, Tissue kallikrein, Kik-9, Klks3
Accession No.	Uniprot ID:P07647
Uniprot	P07647
GeneID	292868;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Kik9

Note: This product is for in vitro research use only