

## Protein lin-28 homolog A

Catalog No: #AP84807



Package Size: #AP84807-1 50ug #AP84807-2 100ug #AP84807-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein lin-28 homolog A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lin-28A, Testis-expressed protein 17, Lin28, Tex17
Accession No.	Uniprot ID:Q8K3Y3
Uniprot	Q8K3Y3
GeneID	83557;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Lin28a

Note: This product is for in vitro research use only