

Lysophosphatidic acid receptor 4

Catalog No: #AP84808



Package Size: #AP84808-1 50ug #AP84808-2 100ug #AP84808-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysophosphatidic acid receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LPA receptor 4,G-protein coupled receptor 23,P2Y purinoceptor 9,P2Y9,P2Y5-like receptor,Purinergic receptor 9,GPR23,LPA4,P2RY9
Accession No.	Uniprot ID:Q99677
Uniprot	Q99677
GeneID	2846;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LPAR4

Note: This product is for in vitro research use only