

Lipid phosphate phosphatase-related protein type 4

Catalog No: #AP84810



Package Size: #AP84810-1 50ug #AP84810-2 100ug #AP84810-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lipid phosphate phosphatase-related protein type 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain-specific phosphatidic acid phosphatase-like protein 1,Plasticity-related gene 1 protein,PRG-1,KIAA0455,PRG1
Accession No.	Uniprot ID:Q7Z2D5
Uniprot	Q7Z2D5
GeneID	9890;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LPPR4

Note: This product is for in vitro research use only