

L-selectin

Catalog No: #AP84813

Package Size: #AP84813-1 50ug #AP84813-2 100ug #AP84813-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD62 antigen-like family member L,Leukocyte adhesion molecule 1,LAM-1,Leukocyte surface antigen Leu-8,Leukocyte-endothelial cell adhesion molecule 1,LECAM1,Lymph node homing receptor,TQ1,gp90-MEL,CD62L,LNHR,LYAM1
Accession No.	Uniprot ID:P14151
Uniprot	P14151
GeneID	6402;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SELL

Note: This product is for in vitro research use only