

cAMP-dependent protein kinase catalytic subunit alpha

Catalog No: #AP84815

Package Size: #AP84815-1 50ug #AP84815-2 100ug #AP84815-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	cAMP-dependent protein kinase catalytic subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PKA C-alpha
Accession No.	Uniprot ID:Q8MJ44
Uniprot	Q8MJ44
GeneID	403556;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKACA

Note: This product is for in vitro research use only