

Lipopolysaccharide-binding protein

Catalog No: #AP84816



Package Size: #AP84816-1 50ug #AP84816-2 100ug #AP84816-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lipopolysaccharide-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LBP
Accession No.	Uniprot ID:Q2TBI0
Uniprot	Q2TBI0
GeneID	512242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LBP

Note: This product is for in vitro research use only