

Lymphoid enhancer-binding factor 1

Catalog No: #AP84817



Package Size: #AP84817-1 50ug #AP84817-2 100ug #AP84817-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphoid enhancer-binding factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LEF-1
Accession No.	Uniprot ID:Q9QXN1
Uniprot	Q9QXN1
GeneID	161452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lef1

Note: This product is for in vitro research use only