

L-lactate dehydrogenase B chain

Catalog No: #AP84820

Package Size: #AP84820-1 50ug #AP84820-2 100ug #AP84820-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-lactate dehydrogenase B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LDH-B,LDH heart subunit,LDH-H,Renal carcinoma antigen NY-REN-46
Accession No.	Uniprot ID:P07195
Uniprot	P07195
GeneID	3945;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LDHB

Note: This product is for in vitro research use only