

Lymphocyte antigen 6 complex locus protein G6f

Catalog No: #AP84821



Package Size: #AP84821-1 50ug #AP84821-2 100ug #AP84821-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphocyte antigen 6 complex locus protein G6f
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C6orf21,G6F,LY6G6D,NG32
Accession No.	Uniprot ID:Q5SQ64
Uniprot	Q5SQ64
GeneID	259215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LY6G6F

Note: This product is for in vitro research use only