

Galectin-3

Catalog No: #AP84826



Package Size: #AP84826-1 50ug #AP84826-2 100ug #AP84826-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gal-3,35 kDa lectin,Carbohydrate-binding protein 35,CBP 35,Galactose-specific lectin 3,IgE-binding protein,L-34 galactoside-binding lectin,Laminin-binding protein,Lectin L-29,Mac-2 antigen
Accession No.	Uniprot ID:P16110
Uniprot	P16110
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lgals3

Note: This product is for in vitro research use only