

Myristoylated alanine-rich C-kinase substrate

Catalog No: #AP84827



Package Size: #AP84827-1 50ug #AP84827-2 100ug #AP84827-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Myristoylated alanine-rich C-kinase substrate |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | MARCKS,Mac |
| Accession No. | Uniprot ID:P26645 |
| Uniprot | P26645 |
| GeneID | 17118; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Marcks

Note: This product is for in vitro research use only