

Mast cell protease 8

Catalog No: #AP84830



Package Size: #AP84830-1 50ug #AP84830-2 100ug #AP84830-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mast cell protease 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rMCP-8, Mast cell protease VIII, rMCP-VIII
Accession No.	Uniprot ID:P97594
Uniprot	P97594
GeneID	29269;54269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mcpt8

Note: This product is for in vitro research use only