

# Integrin alpha-M

Catalog No: #AP84836



Package Size: #AP84836-1 50ug #AP84836-2 100ug #AP84836-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin alpha-M
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD11 antigen-like family member B,CR-3 alpha chain,Cell surface glycoprotein MAC-1 subunit alpha,Leukocyte adhesion receptor MO1,CD11b
Accession No.	Uniprot ID:P05555
Uniprot	P05555
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Itgam

Note: This product is for in vitro research use only