

Integrin beta-1

Catalog No: #AP84839



Package Size: #AP84839-1 50ug #AP84839-2 100ug #AP84839-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CSAT antigen,JG22 antigen,RGD-receptor
Accession No.	Uniprot ID:P07228
Uniprot	P07228
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGB1

Note: This product is for in vitro research use only