

Low-density lipoprotein receptor-related protein 4

Catalog No: #AP84841



Package Size: #AP84841-1 50ug #AP84841-2 100ug #AP84841-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low-density lipoprotein receptor-related protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LRP-4, Multiple epidermal growth factor-like domains 7, KIAA0816, LRP10, MEGF7
Accession No.	Uniprot ID:O75096
Uniprot	O75096
GeneID	4038;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LRP4

Note: This product is for in vitro research use only