

Ligand-dependent nuclear receptor corepressor-like protein

Catalog No: #AP84842

Package Size: #AP84842-1 50ug #AP84842-2 100ug #AP84842-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ligand-dependent nuclear receptor corepressor-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LCoR-like protein
Accession No.	Uniprot ID:Q8N3X6
Uniprot	Q8N3X6
GeneID	254251;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LCORL

Note: This product is for in vitro research use only