

Thymidine kinase 2, mitochondrial

Catalog No: #AP84847



Package Size: #AP84847-1 50ug #AP84847-2 100ug #AP84847-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thymidine kinase 2, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mt-TK
Accession No.	Uniprot ID:O00142
Uniprot	O00142
GeneID	7084;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TK2

Note: This product is for in vitro research use only