

Potassium voltage-gated channel subfamily B member 2

Catalog No: #AP84849

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84849-1 50ug #AP84849-2 100ug #AP84849-3 1mg

Description

Product Name	Potassium voltage-gated channel subfamily B member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Voltage-gated potassium channel subunit Kv2.2
Accession No.	Uniprot ID:A6H8H5
Uniprot	A6H8H5
GeneID	98741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Kcnb2

Note: This product is for in vitro research use only