

Laminin subunit beta-3

Catalog No: #AP84850



Package Size: #AP84850-1 50ug #AP84850-2 100ug #AP84850-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Laminin subunit beta-3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Epiligrin subunit bata,Kalinin B1 chain,Kalinin subunit beta,Laminin-5 subunit beta,Nicein subunit beta |
| Accession No. | Uniprot ID:Q61087 |
| Uniprot | Q61087 |
| GeneID | 16780; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamb3

Note: This product is for in vitro research use only