

Lysosome-associated membrane glycoprotein 1

Catalog No: #AP84854

Package Size: #AP84854-1 50ug #AP84854-2 100ug #AP84854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysosome-associated membrane glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	LAMP-1,Lysosome membrane glycoprotein LEP100
Accession No.	Uniprot ID:P05300
Uniprot	P05300
GeneID	396220;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAMP1

Note: This product is for in vitro research use only