

## Linker for activation of T-cells family member 1

Catalog No: #AP84856



Package Size: #AP84856-1 50ug #AP84856-2 100ug #AP84856-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Linker for activation of T-cells family member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	36 kDa phospho-tyrosine adapter protein,pp36,p36-38
Accession No.	Uniprot ID:O54957
Uniprot	O54957
GeneID	16797;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lat

Note: This product is for in vitro research use only