

Heat shock protein beta-6

Catalog No: #AP84863



Package Size: #AP84863-1 50ug #AP84863-2 100ug #AP84863-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HspB6, Heat shock 20 kDa-like protein p20
Accession No.	Uniprot ID:P97541
Uniprot	P97541
GeneID	192245;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hspb6

Note: This product is for in vitro research use only