

Glucokinase

Catalog No: #AP84864



Package Size: #AP84864-1 50ug #AP84864-2 100ug #AP84864-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glucokinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hexokinase type IV, HK IV, Hexokinase-4, HK4, Hexokinase-D, Gk
Accession No.	Uniprot ID:P52792
Uniprot	P52792
GeneID	103988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gck

Note: This product is for in vitro research use only