

Interferon-stimulated gene 20 kDa protein

Catalog No: #AP84867



Package Size: #AP84867-1 50ug #AP84867-2 100ug #AP84867-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon-stimulated gene 20 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Accession No.	Uniprot ID:Q66UW5
Uniprot	Q66UW5
GeneID	448855;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:lsg20

Note: This product is for in vitro research use only