

## Kielin/chordin-like protein

Catalog No: #AP84869



Package Size: #AP84869-1 50ug #AP84869-2 100ug #AP84869-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kielin/chordin-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cysteine-rich BMP regulator 2,Cysteine-rich motor neuron 2 protein,CRIM-2,Kielin/chordin-like protein 1,KCP-1,CRIM2,KCP1
Accession No.	Uniprot ID:Q6ZWJ8
Uniprot	Q6ZWJ8
GeneID	375616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:KCP

Note: This product is for in vitro research use only