

Kelch-like ECH-associated protein 1

Catalog No: #AP84871



Package Size: #AP84871-1 50ug #AP84871-2 100ug #AP84871-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kelch-like ECH-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytosolic inhibitor of Nrf2, INrf2, Kelch-like protein 19, INRF2, KIAA0132, KLHL19
Accession No.	Uniprot ID:Q14145
Uniprot	Q14145
GeneID	9817;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KEAP1

Note: This product is for in vitro research use only