

# Mediator of RNA polymerase II transcription subunit 31

Catalog No: #AP84875

Package Size: #AP84875-1 50ug #AP84875-2 100ug #AP84875-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mediator of RNA polymerase II transcription subunit 31         |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                                   |
| Species Reactivity | Ms   |
| Other Names        | Mediator complex subunit 31,Mediator complex subunit SOH1,Soh1 |
| Accession No.      | Uniprot ID:Q9CXU1  |
| Uniprot            | Q9CXU1   |
| GeneID             | 67279;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)         |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Med31

Note: This product is for in vitro research use only