

# Integrin alpha-1

Catalog No: #AP84882



Package Size: #AP84882-1 50ug #AP84882-2 100ug #AP84882-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD49 antigen-like family member A,Laminin and collagen receptor,VLA-1,CD49a
Accession No.	Uniprot ID:P56199
Uniprot	P56199
GeneID	3672;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ITGA1

Note: This product is for in vitro research use only