

Kinesin light chain 1

Catalog No: #AP84889



Package Size: #AP84889-1 50ug #AP84889-2 100ug #AP84889-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kinesin light chain 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KLC 1,KLC,KNS2
Accession No.	Uniprot ID:Q07866
Uniprot	Q07866
GeneID	3831;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KLC1

Note: This product is for in vitro research use only