

Leukosialin

Catalog No: #AP84891



Package Size: #AP84891-1 50ug #AP84891-2 100ug #AP84891-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leukosialin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Galactoglycoprotein,GALGP,Leukocyte sialoglycoprotein,Sialophorin,CD43,CD43
Accession No.	Uniprot ID:P16150
Uniprot	P16150
GeneID	6693;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SPN

Note: This product is for in vitro research use only