

Bifunctional arginine demethylase and lysyl-hydroxylase JMJD6

Catalog No: #AP84897

Orders: order@signalwayantibody.com

Package Size: #AP84897-1 50ug #AP84897-2 100ug #AP84897-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Bifunctional arginine demethylase and lysyl-hydroxylase JMJD6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Histone arginine demethylase JMJD6, JmjC domain-containing protein 6, Jumonji domain-containing protein 6, Lysyl-hydroxylase JMJD6, Peptide-lysine 5-dioxygenase JMJD6, Phosphatidylserine receptor, Protein PTDSR, Kiaa0585, Ptdsr
Accession No.	Uniprot ID:Q9ERI5
Uniprot	Q9ERI5
GeneID	107817;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Jmjd6

Note: This product is for in vitro research use only