

Mitogen-activated protein kinase kinase kinase 7

Catalog No: #AP84904



Package Size: #AP84904-1 50ug #AP84904-2 100ug #AP84904-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase kinase kinase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:P0C8E4
Uniprot	P0C8E4
GeneID	100910771;313121;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map3k7

Note: This product is for in vitro research use only