

Lutropin subunit beta

Catalog No: #AP84906



Package Size: #AP84906-1 50ug #AP84906-2 100ug #AP84906-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Lutropin subunit beta |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Lutropin beta chain,Luteinizing hormone subunit beta,LH-B |
| Accession No. | Uniprot ID:P01230 |
| Uniprot | P01230 |
| GeneID | 25329; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lhb

Note: This product is for in vitro research use only