

Arachidonate 12-lipoxygenase, 12S-type

Catalog No: #AP84908



Package Size: #AP84908-1 50ug #AP84908-2 100ug #AP84908-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Arachidonate 12-lipoxygenase, 12S-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	12S-LOX,12S-lipoxygenase
Accession No.	Uniprot ID:P27479
Uniprot	P27479
GeneID	282139;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALOX12

Note: This product is for in vitro research use only