

Tonin

Catalog No: #AP84914



Package Size: #AP84914-1 50ug #AP84914-2 100ug #AP84914-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tonin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Esterase 1, Glandular kallikrein-2, rGK-2, RSKG-5, S2 kallikrein, Klk-2, Ton
Accession No.	Uniprot ID:P00759
Uniprot	P00759
GeneID	24841;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Klk2

Note: This product is for in vitro research use only