

Pyruvate kinase PKM

Catalog No: #AP84916



Package Size: #AP84916-1 50ug #AP84916-2 100ug #AP84916-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyruvate kinase PKM
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pyruvate kinase muscle isozyme, Pk3, Pkm2, Pykm
Accession No.	Uniprot ID:P52480
Uniprot	P52480
GeneID	18746;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pkm

Note: This product is for in vitro research use only