

Epididymal-specific lipocalin-10

Catalog No: #AP84920



Package Size: #AP84920-1 50ug #AP84920-2 100ug #AP84920-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Epididymal-specific lipocalin-10 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Accession No. | Uniprot ID:Q6JVE6 |
| Uniprot | Q6JVE6 |
| GeneID | 414332; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LCN10

Note: This product is for in vitro research use only