

L-lactate dehydrogenase C chain

Catalog No: #AP84921



Package Size: #AP84921-1 50ug #AP84921-2 100ug #AP84921-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | L-lactate dehydrogenase C chain |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | LDH-C,LDH testis subunit,LDH-X,Ldh-3,Ldh3 |
| Accession No. | Uniprot ID:P00342 |
| Uniprot | P00342 |
| GeneID | 16833; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ldhc

Note: This product is for in vitro research use only