

L-lactate dehydrogenase C chain

Catalog No: #AP84922

Package Size: #AP84922-1 50ug #AP84922-2 100ug #AP84922-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-lactate dehydrogenase C chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	LDH-C,LDH testis subunit,LDH-X
Accession No.	Uniprot ID:Q9TSX5
Uniprot	Q9TSX5
GeneID	100502559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LDHC

Note: This product is for in vitro research use only