

## Ly6/PLAUR domain-containing protein 3

Catalog No: #AP84924



Package Size: #AP84924-1 50ug #AP84924-2 100ug #AP84924-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ly6/PLAUR domain-containing protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GPI-anchored metastasis-associated protein C4.4A homolog,C4.4a
Accession No.	Uniprot ID:Q91YK8
Uniprot	Q91YK8
GeneID	72434;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lypd3

Note: This product is for in vitro research use only